

Content Colifast field kit

Rev 150405

- 1 Colifast Micro Detector with storage pouch
- 1 incubator with thermometer. 220V, 44°C
- 1 car converter 12V – 220V

- 18 prefilled vials with Colifast 6 medium (tests)
- 2 vials with 0,5M NaOH solution (developer)
- 1 vial with "blank" (calibration)
- 1 vial with "cal 48" (calibration)
- 1 vial with 10% Na-thiosulfate solution (chloride neutralization)

- 100 cuvettes (disposable)
- 100 cuvette caps (disposable)
- 1 cuvette rack
- 5 small pastettes for addition of NaOH (100µl = 3 drops, can be reused)
- 20 large pastettes for sub-sampling (3ml graded, one for each test)
- 3 sterile 10 ml syringes for addition of water sample (can be reused after wash)
- 4 bags virkon. For disinfection.
- 2 sterile filters. For control samples, fits on syringe

- 1 quick guide
- 10 result sheets
- Extra info



COLIFAST®